Session 2: Working environment

Introduction to Linux for Bioinformatics

Overview of Linux for Bioinformatics

What is the command line


Linux Directory Structure

Basic CLI commands

Packages and Environments Management

HTS File Format

HTS Data Source



Overview of Linux for Bioinformatics

Linux: An open-source (free) operating system widely used in bioinformatics due to its stability, security, and powerful command-line interface.



Why Linux in Bioinformatics?

- Tool availability: Most bioinformatics tools are developed and optimized for Linux.
- Command line tools and scripting: Linux's command-line interface and scripting languages (e.g., Bash) allow for the automation of repetitive tasks, streamlining workflows.
- **Efficiency:** Handles large datasets and complex computations effectively.

How to access a Linux-based system?

Local Desktop









From Windows

Remote connection: SSH Secure Shell protocol







What is the command line?

Command Line Interface* (CLI) is a text-based interface that allows users to interact with the computer operating system (OS) using the keyboard.


A *shell* refers to a program that is an intermediary between the user and the operating system (e.g. Bash, zsh, etc.).

How the CLI shell works?

- **Command Prompt:** The CLI presents a prompt (symbol like \$, >, or %) where the user types commands.
- Commands: User's specific commands that the system interprets and executes such as file management (e.g., copying, moving, deleting files), running programs and scripts.

```
Argument
           Command.
mgem@mgem-VM:~$ mkdir project
mgem@mgem-VM:~$ cp -r project project2
```

- **Command execution:** The shell searches the system's PATH to locate and execute the command file associated with the user's input. The operating system then performs the requested actions.
- Output: The system may return output/results directly in the terminal, such as informational messages, requested data, results of the operation, or an error message if something went wrong.



Command Line Interface (CLI)




Graphical User Interface (GUI)



* Also known as shell, terminal, console, command prompts ... etc.

Linux Directory Structure



```
ıgem@mgem-VM:~$ ls -l /
total 3991648
                                    7 Aug 1 13:29 bin -> usr/bin
lrwxrwxrwx
             1 root root
                                 4096 Feb 26 2024 bin.usr-is-merged
drwxr-xr-x
             2 root root
                                 4096 Aug 19 12:38 boot
             4 root root
drwxr-xr-x
                                 4096 Aug 1 13:35 cdrom
drwxr-xr-x
            2 root root
                                 4120 Sep 11 13:02 dev
drwxr-xr-x 19 root root
drwxr-xr-x 152 root root
                                12288 Aug 19 16:47 etc
                                 4096 Aug 1 13:35 home
drwxr-xr-x
            3 root root
                                    7 Aug 1 13:29 lib -> wsr/lib
lrwxrwxrwx
            1 root root
                                    9 Aug 1 13:29 lib64 -> usr/lib64
lrwxrwxrwx
            1 root root
                                 4096 Apr 8 17:37 lib.usr-is-merged
drwxr-xr-x
            2 root root
                                16384 Aug 1 13:29 lost+found
drwx-----
             2 root root
            1 root vboxsf
                                 4096 Aug 20 11:14 media
drwxrwx---
                                 4096 Jul 21 15:46 mnt
drwxr-xr-x
            2 root root
                                                                  mgem
                                 4096 Aug 19 12:37 opt
drwxr-xr-x
            3 root root
                                                                      db
dr-xr-xr-x 295 root root
                                    0 Sep 11 13:01 proc
                                                                      Desktop
                                 4096 Aug 1 13:36 root
drwx----
            4 root root
                                                                      Documents
           34 root root
                                  960 Sep 11 13:02 run
drwxr-xr-x
                                                                      Downloads
                                    8 Aug 1 13:29 sbin ->
lrwxrwxrwx
            1 root root
                                                                      miniconda3
                                 4096 Mar 31 12:00 sbin.us
drwxr-xr-x
            2 root root
                                                                      Music
                                 4096 Jul 21 15:46 srv
drwxr-xr-x
            2 root root
                                                                      Pictures
            1 root root
                           4087349248 Aug 1 13:29 swapfil
                                                                      Public
           13 root root
                                    0 Sep 11 13:01 sys
dr-xr-xr-x
                                12288 Sep 11 13:07
drwxrwxrwt 16 root root
                                                                      Templates
drwxr-xr-x 12 root root
                                 4096 Jul 21 15:46 usr
                                                                      Videos
drwxr-xr-x 11 root root
                                 4096 Aug 1 13:40 var
```

Absolute path: Complete route to a file or directory from the system root (/home/mgem/Documents) **Relative path:** Location of a file or directory in relation to the current working directory (./Documents)

A single dot (.) refers to the current directory Double dots (..) refers to the directory one level up



Basic CLI commands

Command	Description
whoami	Determine the current username (who am I?)
pwd	Present working directory (where am I?)
ls (-la)	List the content of a directory (with details)
$\operatorname{cd}\operatorname{\textit{pathname}}$	Change directory (folder) in the file system
cd/	Move to the root folder of the file system
cd	Move one level up (one folder) in the file system
ср	Copy a file to another folder
mv	Move a file to another folder
mkdir	Creates a new directory (folder)
cat	Display a file
rm filename	Removes a file
rm -r (or rmdir)	Remove a directory
clear	Clears the CLI window
exit	Closes the CLI window
man command	Shows the manual for a given command

Useful tricks to avoid excessive command typing



- Use copy/paste: ctrl+shift+c (copy) and ctrl+shirt+v (paste).
- Use Up/Down arrow keys: Cycle through recently executed commands.
- Use the TAB key: Autocomplete file/directory name.
- history command: list all recently used commands. Users can copy and paste a desired command to execute it again.
- The wildcard * symbol represents a string of any character of any length.



Default login path: /home/username $\equiv \sim$ e.g., /home/username/Documents $\equiv \sim$ /Documents

Basic CLI commands: File manipulation

> or >> : Redirecting output
| : Piping output

Text files: Human-readable, can be viewed and modified using a text editor

- Text documents (e.g., README files)
- Data in text format (e.g., FASTA, FASTQ, VCF, ...)
- Scripts:
 - Bash scripts (*.sh or *.csh)
 - Python scripts (*.py)
 - R scripts (*.R)

Bash script

A bash script is a file containing a sequence of commands that are executed by the bash program line by line.

It allows you to perform a series of actions, such as navigating to a specific directory, creating a folder, and launching a process using the command line.

Script	t manipu	lation
	o illiallipa.	iauloli

#! /bin/bash	Bash script header
chmod +x script.sh	Give execution permission to a bash script
./script.sh	Run a bash script

-rwxrwxr-x 1 mgem mgem 13 Aug 21 12:42 script.sh



Create and print file	
nano first.txt	Create a text file (.txt) named first
echo "Hello world!" > first.txt	Insert "Hollo word" into <i>first.txt</i> file
cat first.txt	Print content of first.txt file
cat file1.txt file2.txt > file3.txt	Concatenate the contents of the file1/2 into file3.txt
echo "RandomText" >> first.txt	Append the content of <i>first.txt</i> file
grep "o" first.txt	Search for pattern "o" matches
cat first.txt grep "o"	Piping the output of cat into grep command
wc -l first.txt	Count the number of lines in <i>first</i> file
grep "o" first.txt wc -l	Count number of lines with pattern "o"
head first.txt	Display top line of a file
tail first.txt	Display bottom of a file

Compression and Archives	
gzip first.txt	Compress first file into first.txt.gz file
gunzip first.txt.gz	Uncompress first file
zip -r folder.zip folder/	Compress folder to folder.zip file
unzip folder.zip	Uncompress folder folder.zip
tar -cvzf folder.tar.gz folder/	Compress folder into folder.tar.gz
tar -xvzf folder.tar.gz	Extract folder.tar.gz

Packages and Environments Management

In bioinformatics, managing software dependencies and environments is vital for:

- Reproducibility: Ensuring that analyses can be replicated reliably.
- Consistency: Maintaining stable environments across different computational tasks.
- Efficiency: Streamlining workflows and reducing conflicts.

Challenges

- Bioinformatics tools are developed in various programming languages.
- Specific library versions are often required for different tools.
- Manual management of these dependencies can be complex, especially in Linux environments.



Conda is a package manager allows specific versions of programs to be installed, alongside their dependencies. Different sets of programs can be installed in different virtual environments.



Python version X Package 1 Package 2 Package 3 Python version Y Package 1 Package 1 Package 1 Package 2 Package 1 Package 2 Package 2 Package 3 Package 4


Packages and Environments Management




- Simplified package management: Installation, management, and updating of software packages and their dependencies across different programming languages.
- **Environment Management:** Allows to create isolated environments for different projects avoiding conflicts between different software requirements.
- Extensive Package Ecosystem: Provides access to a large collection of pre-built packages for various domains.
- Conda Channels: Serve as repositories for hosting and managing packages. (e.g., Bioconda, conda-forge and general community channels)

Bioconda is a community-enabled repository of 3,000+ bioinformatics packages, installable via the conda package manager.









Bioconda is not available for windows systems

Good news!

Most of the tools required for the analysis pipelines in our hands-on sessions are **pre-installed** within a Conda environment.

The Conda environments can be easily extracted from the virtual machine and transferred to a more powerful computing resource.








Breakdown of the common file formats



Oxford Nanopore Files

What is POD5?

- POD5 is a newer format developed by ONT to store raw signal data from nanopore sequencing runs, replacing the previous FAST5 format
- ONT introduced POD5 to enhance performance in terms of both storage efficiency and access speed
- Contain raw electrical signal data (changes in ionic current) from the nanopore sensor, as well as metadata related to sequencing, like channel and the flow cell Information
- Used by ONT's basecalling software (e.g., Guppy, Dorado), which **translates** the raw signal data into nucleotide sequences (A, T, C, G)



Each base gives a characteristic reduction in the ionic current, allowing the DNA to be sequenced

FASTQ format

- Text-based format used for storing both nucleotide sequences and their corresponding quality score
- Common file extensions: .fastq, and .fq
- Generated directly from sequencing machines (Oxford Nanopore: Fast5 format, PacBio: HDF5 format)

FASTQ file structure

- **1. Header line**: Begins with a @ symbol, followed by a UNIQUE sequence identifier
- 2. Sequence data: The actual nucleotide (A, T, G, C, N) on 1 line
- **3.** A **separator**: Simply as a + sign, sometimes as +SequenceName
- 4. Quality Score: ASCII characters to represent the numerical quality score (Phred +33) for each base Sequence identifier

```
Header line "@" → @d35ee0b2-0e28-44bd-be31-c806b4209a87

Sequence lines → AGGGTTTGATCATGGCTCAGGACGAACGCTGGCGGCGTGCTTAACACATGCAAGTCGAACG

Break → +

Quality score → @=@A?CEDEIFLLJLLMILMMNGCBABCEDACDJMLKJ???>BDDIKKLML@???HDHGH
```

FASTQ format: Understanding Quality String

ASCII TABLE

These characters don't appear in print

Base 33 32 20	e 64
1	
3	nare)
4 4 [END OF TRANSMISSION] 36 24 \$ 68 44 D 100 64 5 5 5 [ENQUIRY] 37 25 % 69 45 E 101 65 6 6 6 [ACKNOWLEDGE] 38 26 & 70 46 F 102 66 7 7 7 [BELL] 39 27 ' 71 47 G 103 67 8 8 8 [BACKSPACE] 40 28 (72 48 H 104 68 9 9 [HORIZONTAL TAB] 41 29) 73 49 I 105 69 10 A [LINE FEED] 42 2A * 74 4A J 106 6A 11 B [VERTICAL TAB] 43 2B + 75 4B K 107 6E 12 C [FORM FEED] 44 2C , 76 4C L 108 6C 13 D [CARRIAGE RETURN] 45 2D - 77 4D M 109 6D 14 E [SHIFT OUT] 46 2E . 78 4E N 110 6E 15 F [SHIFT IN] 47 2F / 79 4F O 111 6F 16 10 [DATA LINE FSCAPE] 48 30 0 80 50 P 112 70 17 11 [DEVICE CONTROL 1] 49 31 1 81 51 Q 113 71 18 12 [DEVICE CONTROL 2] 50 32 2 82 52 R 114 72 19 13 [DEVICE CONTROL 3] 51 33 3 83 53 S 115 73 20 14 [DEVICE CONTROL 4] 52 34 4 84 54 T 116 74 21 15 [NEGATIVE ACKNOWLEDGE] 53 35 5 85 55 U 117 75 22 16 [SYNCHRONOUS IDLE] 54 36 6 86 56 V 118 76 23 17 [END OF TRANS. BLOCK] 55 37 7 87 57 W 119 77 24 18 [CANCEL] 56 38 8 88 88 58 X 120 78 25 19 [END OF MEDIUM] 57 39 9 89 59 Y 121 79	D
5 5 [ENQUIRY] 37 25 % 69 45 E 101 65 6 6 [ACKNOWLEDGE] 38 26 & 70 46 F 102 66 7 7 [BELL] 39 27 ' 71 47 G 103 67 8 8 [BEACKSPACE] 40 28 (72 48 H 104 68 9 9 [HORIZONTAL TAB] 41 29) 73 49 I 105 69 10 A [LINE FEED] 42 2A * 74 4A J 106 6A 11 B [VERTICAL TAB] 43 2B + 75 4B K 107 6E 12 C [FORM FEED] 44 2C , 76 4C L 108 6C 13 D [CARRIAGE RETURN]	c
6 6 [ACKNOWLEDGE] 38 26 & 70 46 F 102 66 7 7 7 [BELL] 39 27 ' 71 47 G 103 67 8 8 8 [BACKSPACE] 40 28 (72 48 H 104 68 9 9 [HORIZONTAL TAB] 41 29) 73 49 I 105 69 10 A [LINE FEED] 42 2A * 74 4A J 106 6A 11 B [VERTICAL TAB] 43 2B + 75 4B K 107 6B 12 C [FORM FEED] 44 2C , 76 4C L 108 6C 13 D [CARRIAGE RETURN] 45 2D - 77 4D M 109 6D 14 E [SHIFT OUT] 46 2E . 78 4E N 110 6E 15 F [SHIFT OUT] 46 2E . 78 4E N 110 6E 16 10 [DATA LINK ESCAPE] 48 30 0 80 50 P 112 70 17 11 [DEVICE CONTROL 1] 49 31 1 81 51 Q 113 71 18 12 [DEVICE CONTROL 2] 50 32 2 82 52 R 114 72 19 13 [DEVICE CONTROL 4] 52 34 4 84 54 T 116 74 21 15 [NEGATIVE ACKNOWLEDGE] 53 35 5 85 55 U 117 75 22 16 [SYNCHRONOUS IDLE] 54 36 6 86 56 V 118 76 23 17 [END OF TRANS. BLOCK] 55 37 7 87 57 W 119 77 24 18 [CANCEL] 56 38 8 88 58 X 120 78 25 19 [END OF MEDIUM] 57 39 9 89 59 Y 121 79	d
7 7 [BELL] 39 27 7 71 47 G 103 67 8 8 [BACKSPACE] 40 28 (72 48 H 104 68 9 9 [HORIZONTAL TAB] 41 29) 73 49 I 105 69 10 A [LINE FEED] 42 2A * 74 4A J 106 6A 11 B [VERTICAL TAB] 43 2B + 75 4B K 107 6E 12 C [FORM FEED] 44 2C , 76 4C L 108 6C 13 D [CARRIAGE RETURN] 45 2D - 77 4D M 109 6E 13 D [CARRIAGE RETURN] 45 2D - 77 4D M 109 6E 15 F [SHIFT IN] 47 2F / 79 4F O 111 6F 16 10 [DATA LINK ESCAPE] 48 30 0 80 50 P 112 70 111 [DEVICE CONTROL 1] 49 31 1 81 51 Q 113 71 18 12 [DEVICE CONTROL 2] 50 32 2 82 52 R 114 72 19 13 [DEVICE CONTROL 2] 50 32 2 82 52 R 114 72 19 13 [DEVICE CONTROL 2] 51 33 3 83 53 S 115 73 20 14 [DEVICE CONTROL 4] 52 34 4 84 54 T 116 74 21 15 [NEGATIVE ACKNOWLEDGE] 53 35 5 85 55 U 117 75 22 16 [SYNCHRONOUS IDLE] 54 36 6 86 56 V 118 76 22 17 [END OF TRANS. BLOCK] 55 37 7 87 57 W 119 77 24 18 [CANCEL] 56 38 8 88 58 X 120 78 25 19 [END OF MEDIUM] 57 39 9 89 59 Y 121 79	e
8 8 [BACKSPACE] 40 28 (72 48 H 104 68 9 9 [HORIZONTAL TAB] 41 29) 73 49 I 105 69 10 A [LINE FEED] 42 2A * 74 4A J 106 6A 11 B [VERTICAL TAB] 43 2B + 75 4B K 107 6E 12 C [FORM FEED] 44 2C , 76 4C L 108 6C 13 D [CARRIAGE RETURN] 45 2D - 77 4D M 109 6E 14 E [SHIFT OUT] 46 2E . 78 4E N 110 6E 15 F [SHIFT IN] 47 2F / 79 4F O 111 6F 16 10 [DATA LINK ESCAPE] 48 30 0 80 50 P 112 70	Ť
9 9 [HORIZONTAL TAB] 41 29) 73 49 I 105 69 10 A [LINE FEED] 42 2A * 74 4A J 106 6A 11 B [VERTICAL TAB] 43 2B + 75 4B K 107 6B 12 C [FORM FEED] 44 2C , 76 4C L 108 6C 13 D [CARRIAGE RETURN] 45 2D - 77 4D M 109 6D 14 E [SHIFT OUT] 46 2E . 78 4E N 110 6E 15 F [SHIFT IN] 47 2F / 79 4F O 111 6F 16 10 [DATA LINK ESCAPE] 48 30 0 80 50 P 112 70 17 11 [DEVICE CONTROL 1] 49 31 1 81 51 Q 113 71 18 12 [DEVICE CONTROL 2] 50 32 2 82 52 R 114 72 19 13 [DEVICE CONTROL 4] 52 34 4 84 54 T 116 74 21 15 [NEGATIVE ACKNOWLEDGE] 53 35 5 85 55 U 117 75 22 16 [SYNCHRONOUS IDLE] 54 36 6 86 56 V 118 76 23 17 [END OF TRANS. BLOCK] 55 37 7 87 57 W 119 77 24 18 [CANCEL] 56 38 8 88 58 X 120 78 25 19 [END OF MEDIUM] 57 39 9 89 59 Y 121 79	g
10	h
11 B [VERTICAL TAB] 43 2B + 75 4B K 107 6E 12 C [FORM FEED] 44 2C , 76 4C L 108 6C 13 D [CARRIAGE RETURN] 45 2D - 77 4D M 109 6C 14 E [SHIFT OUT] 46 2E . 78 4E N 110 6E 15 F [SHIFT IN] 47 2F / 79 4F O 111 6F 16 10 [DATA LINK ESCAPE] 48 30 O 80 50 P 112 70 17 11 [DEVICE CONTROL 1] 49 31 1 81 51 Q 113 71 18 12 [DEVICE CONTROL 2] 50 32 2 82 52 R 114 72 19 13 [DEVICE CONTROL 3] 51 33 3 83 53 S 115	
12 C [FORM FEED] 44 2C , 76 4C L 108 6C 13 D [CARRIAGE RETURN] 45 2D - 77 4D M 109 6C 14 E [SHIFT OUT] 46 2E . 78 4E N 110 6E 15 F [SHIFT IN] 47 2F / 79 4F O 111 6E 16 10 [DATA LINK ESCAPE] 48 30 0 80 50 P 112 70 17 11 [DEVICE CONTROL 1] 49 31 1 81 51 Q 113 71 18 12 [DEVICE CONTROL 2] 50 32 2 82 52 R 114 72 19 13 [DEVICE CONTROL 3] 51 33 3 83 53 S 115 73 20 14 [DEVICE CONTROL 4] 52 34 4 84 54 T 116	Į.
13 D [CARRIAGE RETURN] 45 2D - 77 4D M 109 6D 14 E [SHIFT OUT] 46 2E . 78 4E N 110 6E 15 F [SHIFT IN] 47 2F / 79 4F O 111 6E 16 10 [DATA LINK ESCAPE] 48 30 O 80 50 P 112 70 17 11 [DEVICE CONTROL 1] 49 31 1 81 51 Q 113 71 18 12 [DEVICE CONTROL 2] 50 32 2 82 52 R 114 72 19 13 [DEVICE CONTROL 3] 51 33 3 83 53 S 115 73 20 14 [DEVICE CONTROL 4] 52 34 4 84 54 T 116 74 21 15 [NEGATIVE ACKNOWLEDGE] 53 35 5 85 55 U <t< th=""><th>k</th></t<>	k
14 E [SHIFT OUT] 46 2E . 78 4E N 110 6E 15 F [SHIFT IN] 47 2F / 79 4F O 111 6E 16 10 [DATA LINK ESCAPE] 48 30 O 80 50 P 112 70 17 11 [DEVICE CONTROL 1] 49 31 1 81 51 Q 113 71 18 12 [DEVICE CONTROL 2] 50 32 2 82 52 R 114 72 19 13 [DEVICE CONTROL 3] 51 33 3 83 53 S 115 73 20 14 [DEVICE CONTROL 4] 52 34 4 84 54 T 116 74 21 15 [NEGATIVE ACKNOWLEDGE] 53 35 5 85 55 U 117 75 22 16 [SYNCHRONOUS IDLE] 54 36 6 86 56 V	
15 F [SHIFT IN] 47 2F / 79 4F 0 111 6F 16 10 [DATA LINK ESCAPE] 48 30 0 80 50 P 112 70 17 11 [DEVICE CONTROL 1] 49 31 1 81 51 Q 113 71 18 12 [DEVICE CONTROL 2] 50 32 2 82 52 R 114 72 19 13 [DEVICE CONTROL 3] 51 33 3 83 53 S 115 73 20 14 [DEVICE CONTROL 4] 52 34 4 84 54 T 116 74 21 15 [NEGATIVE ACKNOWLEDGE] 53 35 5 85 55 U 117 75 22 16 [SYNCHRONOUS IDLE] 54 36 6 86 56 V 118 76 23 17 [END OF TRANS. BLOCK] 55 37 7 87 57 W<	m
16 10 [DATA LINK ESCAPE] 48 30 0 80 50 P 112 70 17 11 [DEVICE CONTROL 1] 49 31 1 81 51 Q 113 71 18 12 [DEVICE CONTROL 2] 50 32 2 82 52 R 114 72 19 13 [DEVICE CONTROL 3] 51 33 3 83 53 S 115 73 20 14 [DEVICE CONTROL 4] 52 34 4 84 54 T 116 74 21 15 [NEGATIVE ACKNOWLEDGE] 53 35 5 85 55 U 117 75 22 16 [SYNCHRONOUS IDLE] 54 36 6 86 56 V 118 76 23 17 [END OF TRANS. BLOCK] 55 37 7 87 57 W 119 77 24 18 [CANCEL] 56 38 8 88 58 X </th <th>n</th>	n
17 11 [DEVICE CONTROL 1] 49 31 1 81 51 Q 113 71 18 12 [DEVICE CONTROL 2] 50 32 2 82 52 R 114 72 19 13 [DEVICE CONTROL 3] 51 33 3 83 53 5 115 73 20 14 [DEVICE CONTROL 4] 52 34 4 84 54 T 116 74 21 15 [NEGATIVE ACKNOWLEDGE] 53 35 5 85 55 U 117 75 22 16 [SYNCHRONOUS IDLE] 54 36 6 86 56 V 118 76 23 17 [END OF TRANS. BLOCK] 55 37 7 87 57 W 119 77 24 18 [CANCEL] 56 38 8 88 58 X 120 78 25 19 [END OF MEDIUM] 57 39 9 89 59 Y	0
18 12 [DEVICE CONTROL 2] 50 32 2 82 52 R 114 72 19 13 [DEVICE CONTROL 3] 51 33 3 83 53 S 115 73 20 14 [DEVICE CONTROL 4] 52 34 4 84 54 T 116 74 21 15 [NEGATIVE ACKNOWLEDGE] 53 35 5 85 55 U 117 75 22 16 [SYNCHRONOUS IDLE] 54 36 6 86 56 V 118 76 23 17 [END OF TRANS. BLOCK] 55 37 7 87 57 W 119 77 24 18 [CANCEL] 56 38 8 88 58 X 120 78 25 19 [END OF MEDIUM] 57 39 9 89 59 Y 121 79	p
19 13 [DEVICE CONTROL 3] 51 33 3 83 53 5 115 73 20 14 [DEVICE CONTROL 4] 52 34 4 84 54 T 116 74 21 15 [NEGATIVE ACKNOWLEDGE] 53 35 5 85 55 U 117 75 22 16 [SYNCHRONOUS IDLE] 54 36 6 86 56 V 118 76 23 17 [END OF TRANS. BLOCK] 55 37 7 87 57 W 119 77 24 18 [CANCEL] 56 38 8 88 58 X 120 78 25 19 [END OF MEDIUM] 57 39 9 89 59 Y 121 79	q
20 14 [DEVICE CONTROL 4] 52 34 4 84 54 T 116 74 21 15 [NEGATIVE ACKNOWLEDGE] 53 35 5 85 55 U 117 75 22 16 [SYNCHRONOUS IDLE] 54 36 6 86 56 V 118 76 23 17 [END OF TRANS. BLOCK] 55 37 7 87 57 W 119 77 24 18 [CANCEL] 56 38 8 88 58 X 120 78 25 19 [END OF MEDIUM] 57 39 9 89 59 Y 121 79	r
21 15 [NEGATIVE ACKNOWLEDGE] 53 35 5 85 55 U 117 75 22 16 [SYNCHRONOUS IDLE] 54 36 6 86 56 V 118 76 23 17 [END OF TRANS. BLOCK] 55 37 7 87 57 W 119 77 24 18 [CANCEL] 56 38 8 88 58 X 120 78 25 19 [END OF MEDIUM] 57 39 9 89 59 Y 121 79	t
22 16 [SYNCHRONOUS IDLE] 54 36 6 86 56 V 118 76 23 17 [END OF TRANS. BLOCK] 55 37 7 87 57 W 119 77 24 18 [CANCEL] 56 38 8 88 58 X 120 78 25 19 [END OF MEDIUM] 57 39 9 89 59 Y 121 79	u
23 17 [END OF TRANS. BLOCK] 55 37 7 87 57 W 119 77 24 18 [CANCEL] 56 38 8 88 58 X 120 78 25 19 [END OF MEDIUM] 57 39 9 89 59 Y 121 79	v
24 18 [CANCEL] 56 38 8 88 58 X 120 78 25 19 [END OF MEDIUM] 57 39 9 89 59 Y 121 79	w
25 19 [END OF MEDIUM] 57 39 9 89 59 Y 121 79	×
	ŷ
	y Z
27 1B [ESCAPE] 59 3B ; 91 5B [123 7B	ī
28 1C [FILE SEPARATOR] 60 3C < 92 5C \ 124 7C	ì
29 1D [GROUP SEPARATOR] 61 3D = 93 5D] 125 7D	3
30 1E [RECORD SEPARATOR] 62 3E > 94 5E ^ 126 7E	~
31 1F [UNIT SEPARATOR] 63 3F ? 95 5F 127 7F	[DEL]

FASTQ format: Understanding Quality String

- Quality strings use ASCII values to encode a two-digit integer using a single character
- Keeps files smaller
- Keeps scores in alignment with bases
- Confidence > Phred > Base value adjustment > Number > Character

$$Q = -10 \times \log_{10}(p)$$

p = probability call is not correct

Quality Score (Q)	Probability	Base Call
	(incorrect base call)	Accuracy
0	1/1	0%
10	1 in 10	90%
20	1 in 100	99%
30	1 in 1000	99.90%
40	1 in 10,000	99.99%

Phred	Prob of	Phred+33
score	Error	Ascii
0	1.00000	!
1	0.79433	"
2	0.63096	#
3	0.50119	\$
4	0.39811	%

- The quality line is always the same length as the sequence line
- A single FASTQ file can contain multiple sequences, each with its own header (millions)
- Sequences can be different length

FASTA format

- Text-based format used for representing nucleotide or protein sequences
- Typical file extensions: .fasta, .fa, .fas, .fna, .faa, fsa
- Commonly used for storing and sharing of sequence data
- Compatible with a broad range of bioinformatics tools and software

FASTA file structure

- **1. Header line**: Begins with a > symbol, followed by an identifier and optional description
- 2. Sequence data: The actual nucleotide (A, T, G, C) or amino acid sequence, presented in a plain text format.


Typically, Sequence data can be spread across multiple lengths of length 80, 70, or 60 bases



- A single FASTA file can contain multiple sequences, each with its own header
- Sequences do not need to have the same length

SAM/BAM format

- SAM (Sequence Alignment/Map) files are a standard format used to store aligned sequence data
- BAM (Binary Alignment/Map) is just the binary version of a SAM file, suitable for fast processing
- Hold aligned sequence reads along with metadata, including:
 - The sequence of nucleotides.
 - Quality scores for each base.
 - Alignment data, including reference genome coordinates.
 - Auxiliary tags such as mapping quality, gaps/mismatches, etc.



HTS Data Source



Where to get HTS data?

- Your own Experiment
- NCBI Sequence Read Archive (SRA): A comprehensive repository for raw sequence data from various projects.



(https://www.ncbi.nlm.nih.gov/sra)

• NCBI Genome: A comprehensive repository for genome sequences of various organisms.

(https://www.ncbi.nlm.nih.gov/genome)

• European Nucleotide Archive (ENA): Provides access to raw reads, assembled sequences, and functional information.



(https://www.ebi.ac.uk/ena/browser/home)

• DNA Data Bank of Japan (DDBJ): Offers sequence data submission and retrieval services, including HTS data.



(https://www.ddbj.nig.ac.jp/index-e)

HTS Data Source

How to retrieve HTS Data?

- 1. Find relevant datasets by keywords, organism, study or accession number
- 2. Access the datasets through the database's web interface, or, using command line tools such as sra-tools (NCBI SRA) or ascp (ENA)
- 3. Before downloading, review associated metadata to ensure the dataset's suitability for the analysis (e.g., sample type, sequencing platform).
- 4. Download raw sequence data in FASTQ format using commands such as prefetch, fastq-dump (for SRA)
- 5. Download genome assemblies in FASTA format using commands such as ncbi-genome-download (for NCBI)

Example: ncbi-genome-download -s refseq --formats fasta -A accessions.txt -o output --metadata-table metadata.tsv bacteria

What's Next?

- Prepare for the Hands-On Session
 - Set up the Linux Virtual Machine
 - Ensure all necessary tools and environments are ready
 - Ensure you understand the directory structure where you will be working
 - Get familiar with the command line interface (CLI)

